



Perkin Elmer
Value
Propositions
1121

Process analysis for improved quality



Vanadis[®] cfDNA Platform is a cell-free DNA Analysis Solution Intended for Chromosome Quantification

- The Vanadis[®] cfDNA Platform is composed of 3 instruments:
- The Vanadis Extract[®] instrument for the automated extraction of cell-free DNA from plasma.
- Vanadis Core[®] Instrument for automated cell-free DNA processing.
- Vanadis View[®] Instrument for automated imaging and counting of fluorescent DNA molecules.
- Its system software is for sample tracking, for quality control and for calculating relative amounts of the fluorescent DNA molecules.
- Note: All instruments come with dedicated tables/benches, mandatory for use with the system. Vanadis View instrument additionally comes with System Software.