W.L. Gore
Value
Propositions
1147

Tubing eliminates contamination

Tubing Eliminates Contamination

W.L Gore

- GORE STA-PURE tubing eliminates fluid stream contamination. Low particulation prevents fluid contamination and can extend filter life. Long-term pump testing proves that STA-PURE tubing lasts over 1000 hours at 600rpm or over 6000 hours at 100rpm. It also eliminates the need to rotate tubing through the pump head, saving time and reducing maintenance. The composition also provides flow stability, with virtually no change in flow rate over the life of the tubing.
- GORE CHEM-SURE tubing offers the same benefits as STA-PURE tubing; however, it offers an additional level of chemical resistance. CHEM-SURE tubing